



0.25 Trypsin-EDTA (1X)

For research use only

Catalogue number: BI-1602

Product Description

Trypsin is a serine protease derived from porcine pancreas. It is a single chain polypeptide of 223 amino acid residue with substrate specificity based on positively charged Lysine and Arginine side chains. Trypsin predominantly cleaves peptide chains at the carboxyl sides of Lysine and Arginine, except when either is followed by Proline. It is most commonly used for dissociation and disaggregation of adherent cells. Trypsin can also act to cleave ester and amide linkages of synthetic derivatives of amino acids. Also, ethylenediaminetetraacetic acid (EDTA), a chelating agent is often added to enhance enzymatic activity of trypsin solution. EDTA acts by neutralizing calcium and magnesium ions that enhance cell to cell adhesion. The typical use for this product is in removing adherent cells from culture surfaces.

Specification

This product is 0.25% Trypsin and EDTA in Dulbecco's phosphate buffered saline. Also, it contains sodium bicarbonate and phenol red in HBSS.

Notes

1. Serine protease inhibitors, including DFP, TLCK, APMSF, AEBSEF, and aprotinin will inhibit Trypsin.
2. The concentration of trypsin necessary to dislodge cells from their substrate is dependent primarily on the cell type, the age of the culture and confluency of cell cultures.
3. This product is stored frozen between -10 and -40°C. Repeated cycles of freezing and thawing should be avoided.
4. This product does contain phenol red. Due to shipment on dry ice, there could be significant carbon dioxide buildup in the package. This CO₂ may enter the solution and lower the pH slightly, giving an orange rather than pinkish color. The orange solution will still be suitable for use, or the pH can be adjusted with sodium hydroxide.
5. Incubating cells with too high a trypsin concentration for a long period can damage cell membranes and kill the cells.
6. Solubilizing trypsin or diluting it from a concentrated solution should be done with a buffered salt solution containing no Ca²⁺ or Mg²⁺.
7. Trypsin solutions should be clear of particulates and flocculent material. Do not use if the solution is cloudy or contains a precipitate. Other evidence of deterioration may include degradation of physical or performance characteristics.

How to Use

Before starting the implementation of any protocol, proper concentration of Trypsin-EDTA (0.25% (BI-1602) or 0.05% (BI-1601)) solution and duration of incubation used for dissociation should be determined empirically for individual cell lines. Generally, Time required for dissociation of cells from a surface depends on cell type, cell density, the potency of trypsin, serum concentration in growth medium and time since the last subculture.



1. Remove the spent medium from the culture vessel by aspiration.
 2. Wash the monolayer by adding balanced salt solution without calcium and magnesium (PBS Minus; BI-1401) to the side of the flask opposite the cells.
 3. Rinse the cell sheet by moving the flask gently for 1 to 2 minutes and discard the wash solution.
 4. Add 0.25% Trypsin-EDTA (BI-1602) solution to the side of the flask opposite the cells. The volume should be sufficient enough to completely cover the monolayer of the cells.
 5. Move the flask gently to ensure that the dissociation solution covers the cell sheet.
 6. Incubate the flask at 37°C for 2 to 3 minutes. Monitor the process by observing the flask under an inverted microscope. When dissociation is complete, the cells will be in suspension and appear rounded. In addition to moving gently, flasks of cell lines that are characteristically difficult to remove from substratum may be tapped to expedite removal.
- Note: please note that the exact time needed to dissociate cells will vary according to the cell line. The dissociation process should be monitored closely to avoid cell damage.**
7. Once the cell dissociation is complete add serum containing complete medium to the flask to inhibit the tryptic activity which may further damage the cells.
 8. Disperse the cells into a single cell suspension by pipetting repeatedly.
 9. Count and subculture the cells.

Quality control

- **Appearance:** Red, clear solution
- **pH:** 7.00 -7.60
- **Sterility:** Sterile
- **Shelf life:** 12 months

Storage

- Upon receipt store the product at -20°C in a freezer that is self-defrosting. Once thawed, the product is stable for about 2 weeks at 2-8°C.
- Repeated freezing and thawing reduces enzymatic activity and should be avoided. Once thawed, the solution can be aliquoted in smaller volumes and frozen for future use.

References

1. Treadwell, P. E., and J. D. Ross. "Growth of HeLa cells in human adult and bovine fetal serum medium." *Experimental Biology and Medicine* 111.1 (1962): 197-201.
2. FORM, DAVID M. "Endothelial Cell Proliferation." (1962).

Citations

1. Bafarani, Alireza Hasani, Ziaoddin Mirhoseini, and Farid Heidari. "Achievement to the goat's spermatogonial stem cells from embryos." *Journal of Animal and Poultry Sciences* 3.1 (2014): 08-18.
2. Alsadat, HosseiniAghozbeniElham, ImaniFooladi Abbas Ali, and Nourani Mohammad Reza. "Evaluation of Biocompatible Bioglass/Gelatin Scaffold Enhanced by Mesenchymal Stem Cell for Bone Tissue Regeneration," *Evaluation* 1.1 (2012): 25-30.



3. Aghozbeni, Elham Alsadat Hosseini, Imani Fooladi Abbas Ali, and Mohammad Reza Nourani. "Evaluation of Biocompatible Bioglass/Gelatin Scaffold Enhanced by Mesenchymal Stem Cell for Bone Tissue Regeneration." *Journal of Applied Tissue Engineering* 1.1 (2015).
4. Golafshan, Nasim, Mahshid Kharaziha, and Mohammadhossein Fathi. "Tough and conductive hybrid graphene-PVA: Alginate fibrous scaffolds for engineering neural construct." *Carbon* 111 (2017): 752-763.
5. Golafshan, Nasim, et al. "Nanohybrid hydrogels of laponite: PVA-Alginate as a potential wound healing material." *Carbohydrate polymers* 176 (2017): 392-401.
6. Kamalipooya, Samaneh, et al. "The effects of static magnetic fields on viability and apoptosis in normal and cancerous cells." *Journal of Iranian Clinical Research* 1.3 (2015): 86-90.
7. Gheitanchi, R., M. Kharaziha, and R. Emadi. "Sr-doped forsterite nanopowder: Synthesis and biological properties." *Ceramics International* 43.15 (2017): 12018-12025.
8. Banitaba, Seyedeh Nooshin, et al. "Fabrication of hollow nanofibrous structures using a triple layering method for vascular scaffold applications." *Fibers and Polymers* 18.12 (2017): 2342-2348.
9. Sisakhtnezhad, Sajjad, Mojdeh Heidari, and Ali Bidmeshkipour. "Eugenol enhances proliferation and migration of mouse bone marrow-derived mesenchymal stem cells in vitro." *Environmental toxicology and pharmacology* 57 (2018): 166-174.
10. Nabavinia, Maryam Sadat, et al. "Comparison of Flow Cytometry and ELISA for Screening of Proper Candidate Aptamer in Cell-SELEX Pool." *Applied biochemistry and biotechnology* 184.2 (2018): 444-452.